

## Medtronic's Endurant(TM) II/IIIs Stent Graft System Receives CE Mark for Use with ChEVAR Parallel Graft Technique

December 6, 2016 5:00 AM CT

 Medtronic plc

*Medtronic is First and Only Company to Receive Approval for a ChEVAR Indication, Providing an Off-the-Shelf Option for Patients with Complex Aneurysms*

**DUBLIN – Dec. 6, 2016** – Medtronic plc (NYSE:MDT) today announced that it has received CE (*Conformité Européenne*) Mark for the Endurant™ II/IIIs stent graft system to treat abdominal aortic aneurysm (AAA) patients using a ChEVAR procedure, a parallel graft chimney technique that uses commercially available balloon expandable covered stents combined with a standard aortic stent graft. This expanded indication in the European Union enables the Endurant II/IIIs stent graft system to be used in patients with complex aneurysms with short aortic neck lengths  $\geq 2$  mm, expanded from the prior neck length indication  $\geq 10$  mm.

"Treating aneurysm patients with short aortic necks has been a long-time challenge for clinicians performing endovascular aneurysm repair (EVAR) to treat AAA patients," said Professor Giovanni B. Torsello, M.D., chief of Vascular Surgery, St. Franziskus Hospital, Münster, Germany and co-author of the PROTAGORAS study. "The availability of a standardized approach which increases anatomical applicability will help establish a new standard for patients with complex forms of AAA that may not have been suited for previous procedures."

The CE Mark is supported by a comprehensive review of clinical data from literature using the Endurant II/IIIs stent graft system with the ChEVAR technique. In the flagship PROTAGORAS study, outcomes were tracked with radiologic follow up over a mean of two years. The study used a standardized procedural approach with the Endurant system and balloon expandable covered stents. The results, which were published in the [Journal of Vascular Surgery](#), demonstrated that standardized use of the Endurant II/IIIs stent graft system with ChEVAR in 128 patients is associated with 100 percent technical success, statistically significant aneurysm sac regression ( $p = .001$ ), 95.7 percent primary patency of the chimney grafts and a low incidence of chimney related reinterventions.

"Medtronic is committed to partnering with our clinical community to provide solutions for challenging patients with complex aortic disease," said Daveen Chopra, vice president and general manager of the Aortic business, which is part of the Aortic & Peripheral Vascular division at Medtronic. "The expanded indication for our Endurant II/IIIs stent graft system is a great example of how we can deliver solutions to address unmet clinical needs and improve standard of care for patients who have aneurysms with short neck lengths. We are excited to expand our leadership in EVAR with the first aortic stent graft approved for use with the ChEVAR technique."

The Endurant II/IIIs stent graft system is based on Medtronic's leading Endurant stent graft system, which is selected for nearly one of every two endovascular AAA repairs globally and has resulted in more than 250,000 successful implants. The original Endurant system received the CE (*Conformité Européenne*) Mark in June 2008. The new expanded ChEVAR indication will be initially commercialized in Europe, and then in other countries that recognize the CE Mark approval. In

the U.S., Food and Drug Administration approval for the Endurant Stent Graft System was granted in December 2010. In the U.S., the Endurant II/IIIs stent system is approved for neck lengths  $\geq 10$  mm and  $\leq 60^\circ$  infra-renal angulation and it is not approved for this expanded indication.

In collaboration with leading clinicians, researchers and scientists worldwide, Medtronic offers the broadest range of innovative medical technology for the interventional and surgical treatment of cardiovascular disease and cardiac arrhythmias. The company strives to offer products and services that deliver clinical and economic value to healthcare consumers and providers around the world.

### **About Medtronic**

Medtronic plc ([www.medtronic.com](http://www.medtronic.com)), headquartered in Dublin, Ireland, is among the world's largest medical technology, services and solutions companies – alleviating pain, restoring health and extending life for millions of people around the world. Medtronic employs more than 88,000 people worldwide, serving physicians, hospitals and patients in approximately 160 countries. The company is focused on collaborating with stakeholders around the world to take healthcare Further, Together.

**Any forward-looking statements are subject to risks and uncertainties such as those described in Medtronic's periodic reports on file with the Securities and Exchange Commission. Actual results may differ materially from anticipated results.**

-end-

---

<sup>1</sup> Donas, K. P. et al. The PROTAGORAS study to evaluate the performance of the Endurant stent graft for patients with pararenal pathologic processes treated by the chimney/snorkel endovascular technique. Journal of Vascular Surgery, Volume 63, 1 - 7.

<sup>2</sup> Freedom from chimney graft-related reinterventions was 93.1 percent.